Dissecting the dynamics of protein corona formation on nanoparticles allows reconstructing deep plasma protein concentrations and discovering novel proteoforms

Alexey Stukalov*, Shadi Ferdosi, Yingxiang Huang, Moaraj Hasan, Brittany Lee, Serafim Batzoglou, Asim Siddiqui, and Daniel Hornburg

Blood plasma is a rich and accessible source of known and potential novel **biomarkers**. Even less is known about the **proteoform** composition of these proteins.

known and novel proteoforms.



Nanoparticle Protein Corona









Understanding the NP-protein interactions on the quantitative level opens new possibilities for the nanoparticle-enabled LC-MS/MS analysis.



NP protein corona profiling allows to streamline the design of nanoparticles interrogating specific physicochemical fraction of the proteome occupied by hundreds to thousands of proteoforms.

Reconstructing Protein Abundance Using **NP-Protein Interaction Profiles**

Protein changes in NP coronas accurately follow protein concentration changes between the biosamples, but the dynamic range compression affects absolute abundance estimates.

Neat Abundance



NP-Proteir **Intensity Profile** before

With machine learning applied to the NP-Protein profilles we can reconstruct neat plasma concentrations.

ML Enhances Absolute Neat Protein Intensity Estimates

Using ML to Predict Protein Enrichment in NP

Similar NP-Protein profiles indicate similar Proteograph Assay vs digestion workflow enrichment.



10 15 20 25 log₂(intensity_{neat})

The model predicts that for the neat plasma digestion protocol intensities of the majority of the proteins identified with Proteograph workflow are below the MS detection limit.





NP-protein profiling improves the processing of the retrieved data, including large cohort studies.



Detecting Proteoforms In NP-Peptide Profiles

Alternative splicing, posttranslational cleavage and other modifications are the reasons the same protein can exist in multiple forms – *proteoforms*. These forms may have distinct biological roles and differ in biochemical properties that manifest themselves in different nanoparticle corona intensity profiles of their peptides.



Initial correlation analysis of the peptide-level profiles using the COPF method^{4,5} has identified 61 protein groups composed of several proteoforms, 44 of those groups contained proteoforms not described



- 1. Ferdosi, Stukalov et al, **Advanced Materials** (2022)
- 2. Blume et al, **Nat Comm** (2020)

References

- 3. Ferdosi et al, **PNAS** (2022)
- 4. Donovan et al, **PLoS ONE** (2023)
- 5. Bludau et al, **Nat Comm** (2021)

