


# Deep, unbiased and quantitative mass spectrometry-based plasma proteome analyses of adaptive response to COVID-19 vaccine




**Ting Huang**\*<sup>1</sup>, Alex Rosa Campos<sup>2</sup>, Ramón Díaz Peña<sup>4</sup>, Svetlana Maurya<sup>4</sup>, Khatereh Motamedchaboki<sup>1</sup>, Jian Wang<sup>1</sup>, Laura R. S. Oliveira<sup>3</sup>, Camila I. Alves<sup>3</sup>, Rafael L. Rosa<sup>3</sup>, Yohana P. C. Alves<sup>3</sup>, Lucélia Santi<sup>3</sup>, and Walter O.Beys-da-Silva<sup>3</sup>

<sup>1</sup>Seer, Inc., Redwood City, CA 94065, USA; <sup>2</sup>Plexium, San Diego, CA, USA; <sup>3</sup>Federal University of Rio Grande do Sul, Porto Alegre, Brazil; <sup>4</sup>Sanford Burnham Prebys, La Jolla, CA, USA.

## Proteograph™ workflow coupled with TMTpro 18Plex

COVID-19 vaccines have been extensively used to immunize a large worldwide population during COVID-19 pandemic. However, individuals respond to the vaccine differently, leading to distinct proteome changes in human and vaccination efficiency. Furthermore, our understanding of the underlying molecular mechanisms behind the adaptive response to COVID-19 vaccine remains limited. Here, we used the Proteograph workflow<sup>1</sup> coupled with TMT labeling to quantify the proteome changes in a cohort of 12 volunteers immunized with two doses of COVID-19 vaccine. Our data confirmed that the Proteograph workflow enables identification and quantification of more than 3,000 proteins from human plasma samples in this study. Since 7 volunteers still contracted COVID-19 after two-dose vaccination, our analysis gained further insights into the personalized response to COVID-19 vaccination.




### Methods

Plasma samples were collected before and after first and second vaccination from each volunteer.

Samples were processed by Proteograph™ Product Suite (Seer, Inc.) using five distinct nanoparticles (NPs). Tryptic peptides were then pooled into one single sample for TMT labeling.


A total of 48 samples were allocated into three 18-plex TMTpro mixtures, followed by peptide fractionation by high pH RP and LC-MS/MS analysis, comprised of an EASY nanoLC system coupled to an Orbitrap™ Fusion™ Lumos™ MS equipped with FAIMS Pro Interface (all from Thermo Fisher Scientific). Peptides were separated using a 25-cm analytical C18 Aurora column (IonOpticks) into 24 fractions. A two-hour gradient with three FAIMS Pro compensation voltages was used in this study.

The raw MS spectra data were processed by SpectraMine™ (Biognosys) to generate peptide identification and quantification data. We then used R/Bioconductor package MSstatsTMT for protein quantification, normalization, and differential analysis.




**Figure 1. Proteograph workflow with TMT labeling.** For each biological sample, five peptide mixtures prepared by five distinct nanoparticles of Proteograph™ Assay were pooled into one single sample. Then, 16 pooled samples from 4 volunteers were labeled with TMTpro 18plex and combined to produce a single TMT mixture.


## Proteograph workflow provide valuable insights on personalized human responses to Covid-19 vaccination




**Figure 2. Abundance of identified protein group counts.** Compared with the existing plasma proteomics studies on COVID-19 vaccination<sup>2,3</sup>, the proposed workflow provides much deeper access to human plasma proteome, i.e., 23,372 peptides and 3,094 proteins. Mapping these proteins to the HPPP protein database<sup>4</sup> shows detection across the entire concentration range of the database, 20% of which fell in the low abundance range.





**Figure 3. Hierarchical clustering analysis.** Data shows that samples from the same volunteer were clustered together and demonstrated a high inter-individual variation.



**Figure 4. Differential protein analysis across the cohort.** None of the proteins were differentially expressed after the first dose but 69 proteins became differential after the second dose.



**Figure 5. Individual-specific hierarchical clustering analysis.** Proteins from each volunteer were clustered into significant discrete groups to illustrate various patterns of protein expression changes after vaccination. Most volunteers shared a common N (up-regulation) pattern, where protein abundances were up-regulated by the first dose, decreased after the first dose, and boosted by the second dose.



**Figure 6. Individual-specific pathway enrichment and clustering analysis.** A) Pathway enrichment analysis shows Individuals enriching most gene ontology biological processes (GO-BP) and immune pathways didn't get COVID-19 within 10 months after vaccination. B) Clustering on immune pathways-related proteins separates individuals who didn't get Covid-19.-

Conclusion

Proteograph workflow coupled with TMT labeling quantifies >3,000 protein groups from a cohort of COVID-19 vaccination study with 48 plasma samples.

Covid-19 vaccine related plasma proteomics data show a high between-subject variance.

Individuals with enriched gene oncology and immune pathways didn't get covid within 10 months after vaccination.

References

<sup>1</sup> Blume et al. Nat. Comm. (2020)

<sup>2</sup> Wang et al. Front. Immunol. (2022)

<sup>3</sup> Wang et al. medRxiv (2022)

<sup>4</sup> Schwenk, et al. Journal of Proteome Research. (2017)

