


XT Peptide Quantification Quick Reference





Peptide Quantification

1. Turn on the SP100, [MPE]², chiller unit, and computer

CHECK

- Deck is cleared before loading labware
- Gas supply is ≥ 110 psi
- [MPE]² pressure is set between 105–110 psi

2. Launch  Instrument Control Software (ICS)
3. Prepare consumables from the Proteograph™ XT PQR Labware Kit, and remove the Pierce™ Quantitative Fluorometric Peptide Assay from refrigeration 30 minutes before use
4. On the ICS main menu, select  **Peptide Quantification** and follow the prompts

*See the User Guide for the Proteograph™ XT Assay for full protocol

5. Pipette 100 μ L Peptide Digest Assay Standard to the Peptide Digest Standard tube
6. Pipette 310 μ L Fluorometric Peptide Assay Reagent to each Peptide Assay Reagent tube
7. Add 10 mL Elution Solution into well 1 of the 4-Well Reservoir
8. Add 12 mL Peptide Assay Buffer into well 4 of the 4-Well Reservoir
9. Load remaining labware on the deck, close the front protective cover, and select **Continue** on the ICS

Daily Maintenance

1. From the ICS menu, select  **Maintenance**, then  **Daily Maintenance**

CHECK

- Instrument is turned on
- [MPE]² waste container fill level

2. If needed, empty [MPE]² waste container per lab policy
3. At the **deck and tip waste** stage, select **Yes** to execute



CHECK

- Work deck cleanliness
- Teaching needles are on the work deck

4. If necessary, perform Weekly Maintenance

Daily Maintenance ensures proper 8-channel pipette function

Weekly Maintenance

1. Ensure the instrument is turned off
2. Clean the tip eject sleeves, stop disk, O-rings, and CO-RE head pipette channels with a lint-free wipe dampened with deionized water or microcide
3. Clean the carriers, tip eject plate, and barcode sensor with a lint-free wipe dampened with 70% isopropyl or ethanol (remove carriers and tip eject plate for cleaning)
4. Clean both sides of front protective cover with a lint-free wipe dampened with deionized water
5. Turn on the instrument, then from the ICS menu select  **Maintenance**, then  **Weekly Maintenance**
6. Follow the ICS prompts

